

ProVetScan SV-3

Reproductive control in porcine,
ovine and caprine





Provetscan SV-3 Wi-Fi Probe



Ultrasound scanner for reproductive diagnostic in Porcine, Ovine and Caprine.

Specifications:

Monitor: Smartphones and tablets of Android and iOS.
Frequency: 3,5 MHz
Scanning mode: Sectorial
Number of crystals: One annular piezo crystal
Scanning depth: 100-180 mm
Wi-Fi chip: 2,4 - 5G
Grey levels: 256 niveles
Battery working time: 5 hours
Size: 140x45x40 mm
Weight: 200 g



Included:

Provetscan SV-3 Wi-Fi probe
Tablet Samsung Galaxy Tab A 8"
Protection cover
Neck belt
Transport suitcase

Optional:

Wireless charger

New Ve Tec S.L.